

# BIOP PTS Instructions Manual

---

Installation & Protocol

T. Laroche O. Burri, A. Seitz

v1.1, June 2016

# Summary

---

1. Installation
2. Interface
3. Settings
4. Performing PTS Analysis on a single image
5. Batch PTS Analysis
6. Summarizing Results
7. Extra notes

# Installation

## Prerequisites:

1. ActionBar plugin from Jérôme Mutterer
2. BIOPLib Installer and BIOPLib.ijm

Both of these can be installed via Fiji's "Update Sites":


1. From Fiji, go to *Help > Update...*
2. Click "*Manage Update Sites*"
3. Check both 'IBMP-CNRS' and 'PTBIOP' update sites
4. Restart Fiji.

## Installation:

1. Get the latest BIOP\_PTS.jar at <https://github.com/lacan/ijm-pts>
2. Copy it inside the 'BIOP' folder in your 'Plugins' folder
3. Restart Fiji


## Testing:

Run the PTS by going to *Plugins > BIOP > BIOP PTS*



# Interface

- The interface is provided as an ActionBar.



# Settings

## = Main settings =

### Oversample Image

Rescales the image 2x in XY to allow for smoother detection by the peak detector

### From Radius

Initial radius [in um] for the detection

### Till Radius

Last radius [in um] for the detection

### Increment

Step size [in um] between radii for analysis

### Line Width

When the circle is straightened, interpolate the data around the circle by this much [pixels]

### Distance Allowed

### Variation Percent:

If the calculated Peak-To-Peak distance lies within the theoretical Peak-To-Peak distance plus or minus this percentage, then it is considered as a valid distance

## = Center Tuning =

### Fine-Tune Center

Attempt to correct for user-set center position by using "Find Maxima" with the parameters below

### Dark Center

For V1 of the PTS, specifies whether the center point is dark on a light background.

### Center Noise Tolerance

### User Draws ROI

Noise tolerance for the "Find Maxima Plugin". The user defined the center manually, otherwise it assumed the target is well centered.

### Measure on

### Flattened Image

Performs pseudo-flatfield correction on the image before analysing

## = Summary Options =

When all three criteria below are met, we consider that this is the OQP for the current image.

### Min Detected Peaks

How many peaks should be detected as a minimum to consider we have found the maximum resolution.

### Min Good Peaks Percent



How many among the peaks should have a distance that matches the theoretical distances

### Max Allowed Delta

### Peak-to-Peak

How much variance should we allow in the detected positive good distances

You can find a setting file with sensible defaults to get started with at <https://github.com/lacan/ijm-pts>




### Settings

Saving and loading analysis settings

### Define Parameters

Dialog with all relevant settings for analysis

Clicking 'OK' generates a text window that should be kept open




You can save these settings for use next time by clicking on "Save Parameters"




# Single Image Workflow

After defining the PTS' settings and opening an image.

1. Ensure image calibration is correct
2. Click on **"Current Image"**
3. If *"User Draws ROI"* is set, you will be presented with an interface to define the center and extents.



Run PTS Analysis  
On Image or Folder




## Image Calibration


Define Objective, CCD Pixel size, and c-mount to establish calibration (if needed)

# Single Image Workflow – Manual Extents

1. Hover over the center cross with the mouse and wait until it changes color.
2. Left click and hold to drag the PTS center
3. Do the same for the two circles that represent from where till where the analysis should be performed



Once the extents are set, right-click to begin the analysis



Run PTS Analysis  
On Image or Folder

**Image Calibration**  
Define Objective, CCD Pixel size, and c-mount to establish calibration (if needed)

# Single Image Workflow – Results

A Result Table called "Pattern Detector Detailed Results" is created with all relevant data:

## Label

Image name

## Circle Radius

Radius of the circle where OQP was measured

## Oversampled

Whether the image was oversampled

## Expected Mean Peak-to-Peak Distance

Theoretical P2P distance, calculated as  $R \sin\left(\frac{2\pi}{360}\right)$   
With R: Circle Radius

## Mean Peak-to-Peak Distance

Calculated P2P distance averaged from all "good" detections

## StDev

Standard Deviation of all good P2P distances

## Expected FWHM, FWHM


The full width at half maximum taken as  $2/3 * P2P$

## Total Peaks

Total number of peaks detected by 1D Maxima finder

## Good Distance %

The amount of peaks that had a P2P distance within bounds defined by "Variation Percent"



File	Edit	Font
Label	Circle Radius [microns]	
1	20160606_lsm710_63x1.4_blue_pol 1.czi	12
2	20160606_lsm710_63x1.4_blue_pol 1.czi	12.500
3	20160606_lsm710_63x1.4_blue_pol 1.czi	13
4	20160606_lsm710_63x1.4_blue_pol 1.czi	13.500
5	20160606_lsm710_63x1.4_blue_pol 1.czi	14
6	20160606_lsm710_63x1.4_blue_pol 1.czi	14.500
7	20160606_lsm710_63x1.4_blue_pol 1.czi	15
8	20160606_lsm710_63x1.4_blue_pol 1.czi	15.500
9	20160606_lsm710_63x1.4_blue_pol 1.czi	16
10	20160606_lsm710_63x1.4_blue_pol 1.czi	16.500
11	20160606_lsm710_63x1.4_blue_pol 1.czi	17
12	20160606_lsm710_63x1.4_blue_pol 1.czi	17.500
13	20160606_lsm710_63x1.4_blue_pol 1.czi	18

## Max Consecutive Good Distances

How many consecutive good distances were detected

## Michelson Contrast


Contrast metric defined as  $\frac{I_{\max} - I_{\min}}{I_{\max} + I_{\min}}$   
Where  $I_{\max}$  and  $I_{\min}$  are the averaged minimum and maximum intensities

found by the 1D peak finder



# Checking Analysis

After analysis, you can click on a circle in the ROI Manager and see how the analysis was performed on the bright and dark stripes



**Check Peaks**  
Shows the current ROI  
straightened and displays  
the good peak detections

The color bars indicate what is considered as a "Good Peak to Peak Distance", as defined by the theoretical values plus or minus the **"Variation Percent"**

Anything shorter than **blue** is too short and longer than **red** is too long.

# Batch Processing


By clicking on “**Choose Folder**” you will be prompted for a folder.

After selection, each image will be opened and analysed in the same way as before.

Once an image is finished, it will be closed and the next one will open.

The results will all be appended to the table “**Pattern Detector Detailed Results**”

The ROI sets are saved for each image in a folder called “ROI Sets”




# Summary Table

The current implementation displays all analyzed values. To find the one matching the OQP, one can use **"Summarize Table"**

A new table containing a single line per image is produced.

This is obtained by finding for each image, the first line in the Detailed results matching the *"Summary Options"* defined under **"Detection Settings"**

Get the best OQP per image  
Based on Detailed Results  
Table



Pattern Detector Summary


—

□

✕

File Edit Font

is]	StDev	Expected FWHM [microns]	FWHM [microns]	Total Peaks	Good Distance [%]	M ▲
0.015	0.221		0.221	358	99.721	21
0.016	0.227		0.226	358	100	3:
0.009	0.192		0.194	356	95.506	9:



## Note

To be able to summarize, the table called **"Pattern Detector Detailed Results"** must be open.

# Extra Information From File Name

You can get extra information using proper filenaming convention

Example: *20150330\_Zeiss-IN2\_63x14\_white\_damage 1.lsm*

Would export the following:

**Date** : 20150330

**Machine** : Zeiss-IN2

**Objective** : 63x

**NA** : 1.4

**wavelength** : white

**Note** : damage

Pattern Detector Detailed Results

File Edit Font

iled	Date	Machine	Objective	NA Objective	Wavelength	Note	Expected Me
	20160606	Ism710	63x	1.4	blue	pol	0.401
	20160606	Ism710	63x	1.4	blue	pol	0.410
	20160606	Ism710	63x	1.4	blue	pol	0.419
	20160606	Ism710	63x	1.4	blue	pol	0.428
	20160606	Ism710	63x	1.4	blue	pol	0.436
	20160606	Ism710	63x	1.4	blue	pol	0.445
	20160606	Ism710	63x	1.4	blue	pol	0.454
	20160606	Ism710	63x	1.4	blue	pol	0.462

These columns can help make sense of the data when compiling large amounts of PTS images